




Guide to self-testing for Coronavirus using a




Lateral Flow Device


 We  are taking  this  test  to  keep safe  and  to keep

 people  around  you  safe  too.




 You can  watch a video  about  the test


<https://learninghub.nhs.uk/self-swab>



Before



you



start

1



Wash your hands with soap for 20 seconds.

2



Dry your hands

3



Clean and dry a flat surface, where you will take the test.

4



Check that you have everything you need for one test



Nozzle cap



Extraction buffer bottle

Swab, inside sealed wrapper



Extraction tube



Rapid test strip




5



write down the lot number.



Prepare



your



test

1




Gently blow your nose.

2




Wash and dry your hands again
(or use sanitiser if available)

3



check the swab is sealed, find the
soft, fabric tip of the swab.

4



Place the extraction tube in a
small cup to keep it upright

5



Unscrew the extraction buffer pot

6



Try to avoid touching the tube
with the bottle.

7




Place 6 drops only of the extraction
buffer into the extraction tube.


8



Reseal the extraction buffer for
further use.



Take



your



nose swab

9



Open the package, take out the swab. **DON'T** touch the fabric swab.

10



Put the swab into one nostril and roll around your nostril **10** times slowly

11



Now do the same in the other nostril.

12



Put the fabric end of the swab into the extraction tube. Press the swab on the edge and twist around the tube for 10 seconds.

13



Take out the swab while squeezing the tube and fabric swab.

14



Put the swab back into the packet and put it in the bin.

15



Press the cap tightly onto the tube to avoid any leaks.

16




Turn the tube upside down and lightly squeeze 2 drops into the specimen well

17




Set a timer for 30 minutes. Leave the test strip on a flat surface.

18



Do not read the results until 30 minutes have passed.



Your

results



Positive result

Two lines - even faint lines is positive



Negative result

This shows a negative result



Invalid test

The test has failed and needs to be retaken



Report



your



positive



test



result

